## **SPEIKO® Zink-Phosphat-Zement (vormals Zhanelka Microphos)**

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Pulver: Zinkoxid, Siliciumdioxid, Aluminiumoxid, Calciumoxid, Natriumoxid, Magnesiumoxid.

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Inhalt: 20 ml Art. Nr. 1457

Stand der Information: November 2005

## **SPEIKO® Zinc phosphate cement (previously Zhanelka Microphos)**

### **Composition:**

Powder: zinc oxide, silicon dioxide, aluminium oxide, calcium oxide, sodium oxide, magnesium oxide.  
SPEIKO® Phospho universal liquid for normal-setting and quick-setting: aluminium and zinc in phosphoric acid.

### **Applications:**

Pulp protection material under metal and synthetic fillings, and as a fixing cement for crowns and bridges, dowels and orthodontic bands. Also suitable for use as a temporary filling material. Build-up material.

### **Contraindications:**

As there is a risk of damage to the dental pulp, the product must not be applied to dentine in close proximity to the pulp or to the exposed pulp. Wear must be expected to increase in the course of time.

### **Side effects / interactions:**

Pulp irritation may occur after cementing and with use as pulp protection or as a temporary filling.

### **Dosage instructions, form and duration of application:**

Mix powder in a volumetric ratio of 2.2 : 1 with SPEIKO® Phospho universal liquid, adding the powder to the liquid and not vice versa. The mixing operation should be continued for ca 1 minute until a smooth, pasty consistency is obtained.

Like other dental cements, SPEIKO® zinc phosphate cement is mixed on a glass plate, this should be cooled before preparing larger quantities to prolong the handling time.

All surfaces of the tooth and of the cavity must be prepared, cleaned and dried in advance, so that the freshly mixed cement can be applied. The setting time after application is ca. 2,5 - 8 min., depending on whether normal-setting or quick-setting SPEIKO® Phospho universal liquid is used.

Correct setting can be controlled if required by examining the residues on the glass plate.

### **Special instructions:**

The liquid contains phosphoric acid and has corrosive effect. In the event of inadvertent contact with the skin, mucosa or eyes, rinse immediately with plenty of water. Lack of appropriate pulp protection may result in pulp necrosis. The patient should not eat for the first 2 hours after placement of the filling. This medical product should not be used after the marked date.

Keep out of reach of children.

### **Available in the following packs:**

SPEIKO® zinc phosphate Powder shade 2, light yellow	Contents:	100 g	Art. no. 1455
SPEIKO® Phospho universal liquid, normal-setting	Contents:	20 ml	Art. no. 1454
SPEIKO® Phospho universal liquid, quick-setting	Contents:	20 ml	Art. no. 1457

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**SPEIKO®- Dr. Speier GmbH Postfach 7204 D-48039 Münster**