

DATE OF ISSUE: May 28, 2015

ISSUED BY PASCAL

**Product Name:** Gel Cord<sup>®</sup> / Gel Cord<sup>®</sup> Clear**1. IDENTIFICATION OF THE MATERIAL AND SUPPLIER**

Product Name: Gel Cord<sup>®</sup> / Gel Cord<sup>®</sup> Clear  
Product Use: Topical hemostatic gel  
Manufacturer Name: Pascal International, Inc.  
Address: 2929 NE Northup  
Bellevue, Washington 98004  
Emergency Telephone (Infotrac): 1-800-535-5053, 1-352-323-3500  
Telephone: 1-425-827-4694, 1-800-426-8051  
Fax: 1-425-827-6893

**2. HAZARDS IDENTIFICATION****GHS Classification & Label Elements:****Health Hazard****Mild Eye Irritant** Category 2B**Label Elements:** No symbol**Signal Word:** **Warning****Hazard Statements:**

H320 Causes eye irritation

**Prevention Statements**

P264 Wash thoroughly after handling

P280 Wear protective gloves/ protective clothing/eye protection

**Response Precautionary Statements:**

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing.

**3. COMPOSITION / INFORMATION ON INGREDIENTS**

Chemical Entity (Components)	CAS No.	Proportion
Aluminum Sulfate, USP	10043-01-3	25%
Amorphous Silica	7631-86-9	2.5%

**4. FIRST AID MEASURES**

Inhalation: N/A  
Ingestion: Induce vomiting. Seek medical assistance.  
Skin: Flush well with water for at least 15 minutes.

DATE OF ISSUE: May 28, 2015

ISSUED BY PASCAL

**Product Name:** Gel Cord® / Gel Cord® Clear

Eye: Flush well with water for at least 15 minutes.  
First Aid Facilities: Seek Medical assistance for further treatment, observation and support if necessary.

## 5. FIRE FIGHTING MEASURES

Extinguishing Media: Suitable extinguishing media for the surrounding fire should be used.  
Specific Hazards: None  
Precautions in connection with fire: Use of self contained breathing apparatus is recommended. Wear protective clothing to prevent contact with skin and eyes.

## 6. ACCIDENTAL RELEASE MEASURES

Emergency Procedures: Keep personnel away from area. Contain spills.  
Clean Up Procedure: Clean area with soapy water.

## 7. HANDLING AND STORAGE

Handling: Use original container.  
Storage: Store in cool, well ventilated area.

## 8. EXPOSURE CONTROLS / PERSONAL PROTECTION

Exposure Limits: Not determined  
Respiratory Protection: N/A  
Eye Protection: Safety glasses/goggles  
Hand Protection: Protective gloves  
Engineering controls: Standard room ventilation

## 9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance: Blue gel  
Odor: Raspberry  
pH: 3.0 to 4.0  
Vapor Pressure: Not determined  
Vapor Density: Not determined  
Boiling Point: Not determined  
Freezing / Melting Point: Not determined  
Solubility in Water: Aluminum sulfate is highly soluble in water.  
Specific Gravity: Not determined  
Flash Point: Not determined  
Flammability: N/A

DATE OF ISSUE: May 28, 2015

ISSUED BY PASCAL

**Product Name:** Gel Cord® / Gel Cord® Clear

## 10. STABILITY AND REACTIVITY

Chemical Stability: Stable under ordinary conditions of use and storage  
Conditions to Avoid: High Heat  
Incompatible Materials: None known  
Hazardous Decomposition Products: None known  
Hazardous Reactions: None known

## 11. TOXICOLOGICAL INFORMATION

Hazardous Ingredient: Aluminum Sulfate Mixture  
LD<sub>50</sub> (orally) MICE: 24000 mg/kg  
Routes of Exposure: Refer to Section 4 of SDS

## 12. ECOLOGICAL INFORMATION

Ecotoxicity: No known effects  
Persistence and Degradability: No information found

## 13. DISPOSAL CONSIDERATIONS

Disposal Method and Containers Empty containers may be disposed of in waste containers.  
Special Precautions Disposal should be carried out in compliance with federal, state, and local regulations.

## 14. TRANSPORT INFORMATION

Not a dangerous good. Not regulated  
Hazchem code: None allocated

## 15. REGULATORY INFORMATION

Precautionary Phrases: For professional in-office use only.  
Keep out of reach of children.  
Poison Schedule: None allocated  
AICS (Australia) All ingredients are listed  
Note: Refer to applicable international, national, and local regulations and provisions.

DATE OF ISSUE: May 28, 2015

ISSUED BY PASCAL

<b>Product Name:</b>	<b>Gel Cord® / Gel Cord® Clear</b>
----------------------	------------------------------------

## 16. OTHER INFORMATION

**Legal Disclaimer:** Information contained in this SDS is based on data believed to be correct. This information does not purport to be all inclusive and shall be used only as a guide. No warranty is expressed or implied regarding the accuracy of this data or the results to be obtained from the use thereof. In no event will the manufacturer or the distributor be responsible for damages of any nature whatsoever resulting from the use of or reliance upon this information.

**SDS History**

Date of Issue:	May 28, 2015
Replaces Issue:	June 15, 2014
Changes made to previous issue:	Added GHS classification, pictograms, hazard & precautionary statements